



PeriData Information

Facility _____

completed and checked - date _____ by _____

<p>How do I start PeriData out of my medical software ?</p>	<ul style="list-style-type: none"> ▶ open the patient file <p><i>fill in the steps needed in your medical software:</i></p> <div style="border: 1px dashed gray; height: 100px; margin: 10px 0;"></div> <ul style="list-style-type: none"> ▶ the PeriData window shows all previously stored examinations
<p>How do I send the current examination from the device to PeriData ?</p>	<div style="border: 1px dashed gray; padding: 5px; margin-bottom: 10px;"> <p>start PeriData out of the patient file >> waits for data</p> </div> <ul style="list-style-type: none"> ▶ at the end of the examination, the data are sent automatically when saved on the device <p><i>?? cross this out if it is not automatically and fill in the steps needed in your perimeter device:</i></p> <div style="border: 1px dashed gray; height: 50px; margin-top: 10px;"></div>
<p>How do I send stored examinations from the device to PeriData ?</p>	<ul style="list-style-type: none"> ▶ click on the red perimeter button >> a window opens <div data-bbox="1241 1603 1458 1706" data-label="Image"> </div> <ul style="list-style-type: none"> ▶ click on the button ? HELP >> a detailed description for your device is shown
<p>Where do I call if there are problems ?</p>	<ul style="list-style-type: none"> ▶ +49 22 044 264 770 ▶ email: info@peridata.com